

bs-11937R

[Primary Antibody]

BioSS[™]
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

Protein atonal homolog 8 Rabbit pAb

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 84913</p> <p>Target: Protein atonal homolog 8</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human HATH6/Protein atonal homolog 8: 231-321/321.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: Putative transcription factor. May be implicated in specification and differentiation of neuronal cell lineages in the brain. May participate in kidney development and may be involved in podocyte differentiation.</p>	<p>Isotype: IgG</p> <p>SWISS: Q96SQ7</p>	<p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Pig, Sheep, Cow, Dog)</p> <p>Predicted MW.: 35 kDa</p> <p>Subcellular Location: Nucleus</p>
--	--	--