
Recombinant Human FOXO6 Protein, N-His-SUMO

Catalog Number: bs-105505P

Species: Human

AA Seq: 92-172/492

Predicted MW: 21.65 kDa

Tags: N-His-SUMO

Activity: Not tested

Purity: >90% as determined by SDS-PAGE.

Purification: AC

Form: Lyophilized

Storage: Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.

Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Background: Murine Foxo6 is a member of the murine forkhead family of transcription factors. This family consists of over 30 members, the vast majority of which is important in embryonic development. These forkhead transcription factors may play a role in maintenance and survival of developing and adult neurons.