

---

## Recombinant Human TARSL2 Protein, N-His

Catalog Number: bs-105518P

Species: Human

AA Seq: 158-383/802

Predicted MW: 28.01 kDa

Tags: N-His

Activity: Not tested

Purity: >90% as determined by SDS-PAGE.

Purification: AC

Form: Lyophilized

Storage: Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.

Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Background: TARSL2 belongs to the class-II aminoacyl-tRNA synthetase family and it has two isoforms.