

---

## Recombinant Human MFSD2B Protein, N-GST & C-His

Catalog Number: bs-105521P

Species: Human

AA Seq: 2-112/504

Predicted MW: 39.46 kDa

Tags: N-GST & C-His

Activity: Not tested

Purity: >90% as determined by SDS-PAGE.

Purification: AC

Form: Lyophilized

Storage: Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.

Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Background: MFSD2B (Major Facilitator Superfamily Domain Containing 2B) is a Protein Coding gene. An important paralog of this gene is MFSD2A.