
Recombinant Human MCIDAS Protein, N-GST & C-His

Catalog Number: bs-105526P

Species: Human

AA Seq: 314-385/385

Predicted MW: 36.13 kDa

Tags: N-GST & C-His

Activity: Not tested

Purity: >90% as determined by SDS-PAGE.

Purification: AC

Form: Lyophilized

Storage: Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.

Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Background: Transcription regulator required for multiciliate cell differentiation. Acts by promoting transcription of genes required for multiciliate cell formation. Probably acts in a multiprotein complex (By similarity). Plays a role in mitotic cell cycle progression by promoting cell cycle exit.