

---

## **Recombinant Human CSN3 Protein, N-His**

Catalog Number: bs-105733P

Species: Human

AA Seq: 24-182/182

Predicted MW: 20.12 kDa

Tags: N-His

Purity: >90% as determined by SDS-PAGE.

Purification: AC

Form: Lyophilized

Storage: Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.

Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Background: Important role in determination of the surface properties of the casein micelles.