

---

## Recombinant Human CALCB Protein, N-His

Catalog Number: bs-105742P

Species: Human

AA Seq: 1-127/127

Predicted MW: 16.02 kDa

Tags: N-His

Purity: >90% as determined by SDS-PAGE.

Purification: AC

Form: Lyophilized

Storage: Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.

Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

**Background:** CALCB (Calcitonin Related Polypeptide Beta) is a Protein Coding gene. Diseases associated with CALCB include Immunodeficiency 17 and Thyroid Carcinoma, Familial Medullary. Among its related pathways are Peptide ligand-binding receptors and Presynaptic function of Kainate receptors. Gene Ontology (GO) annotations related to this gene include hormone activity and neuropeptide hormone activity. An important paralog of this gene is CALCA.