

---

## Recombinant Human CD52 Protein, C-Fc

Catalog Number: bs-105772P

Species: Human

AA Seq: 1-61/61

Predicted MW: 30.17 kDa

Tags: C-Fc

Purity: >90% as determined by SDS-PAGE.

Purification: AC

Form: Lyophilized

Storage: Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.

Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Background: Involved in positive regulation of cytosolic calcium ion concentration. Predicted to be located in extracellular region and plasma membrane. Predicted to be intrinsic component of plasma membrane. Predicted to be active in sperm midpiece. [provided by Alliance of Genome Resources, Apr 2022]