
Recombinant Streptococcus pyogenes EndoS Protein, N-GST

Catalog Number: bs-106843P

Species: Streptococcus pyogenes

AA Seq: 37-995/995

Predicted MW: 134.61 kDa

Tags: N-GST

Activity: Not tested

Purity: >90% as determined by SDS-PAGE.

Purification: AC

Form: Lyophilized

Storage: Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.

Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.