

---

## **Recombinant Human CD307e/FCRL5 Protein, C-His**

Catalog Number: bs-107456P

Species: Human

AA Seq: 1-844/977

Predicted MW: 93.24 kDa

Tags: C-His

Activity: Not tested

Purity: >90% as determined by SDS-PAGE.

Purification: AC

Form: Lyophilized

Storage: Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.

Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.