
Human CD40L Mouse mAb

Cat. No:	V7902
Name:	Human CD40L Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	10A11
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Applications:	ELISA
Antigen:	Recombinant Human CD40L Protein
Pairing Information:	
V7902-V7903&V7904	