
Human CD40L Mouse mAb

Cat. No:	V7903
Name:	Human CD40L Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	1H2
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS(pH7.4)
Applications:	ELISA
Antigen:	Recombinant Human CD40L Protein
Pairing Information:	
V7901&V7902-V7903	