
phospho-Tau (Thr181) Mouse mAb, FITC conjugated

Cat. No:	V5227-FITC
Name:	phospho-Tau (Thr181) Mouse mAb, FITC conjugated
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	1B2
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	Human
Concentration:	1mg/ml
Buffer:	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% glycerol.
Applications:	Excitation spectrum: 495nm Emission spectrum: 519nm Not yet tested in other applications. Optimal working dilutions must be determined by the end user.
Antigen:	Recombinant Human p-Tau (Thr181) Protein