
beta-Amyloid (1-40) Mouse mAb, FITC conjugated

Cat. No:	V5312-FITC
Name:	beta-Amyloid (1-40) Mouse mAb, FITC conjugated
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	3A6
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	(predicted: Human)
Concentration:	1mg/ml
Buffer:	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% glycerol.
Antigen:	KLH conjugated synthetic peptide derived from human beta-Amyloid(1-40)