
Human CHI3L1 Recombinant Mouse mAb

Cat. No:	V7801
Name:	Human CHI3L1 Recombinant Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Recombinant
Clone No:	3C5
Host/Source:	Mouse
Subtype:	IgG2a
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS(pH7.4)
Applications:	ELISA
Antigen:	Recombinant Human CHI3L1 Protein (bs-43148P)
Pairing Information:	
V7802&V7803-V7801	