

---

## Human GFAP Recombinant Rabbit mAb

Cat. No:	V5607
Name:	Human GFAP Recombinant Rabbit mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Recombinant
Clone No:	3B7
Host/Source:	Rabbit
Subtype:	IgG
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS(pH7.4)
Applications:	ELISA
Antigen:	Recombinant Human GFAP Protein ( <a href="#">bs-42001P</a> )
Pairing Information:	
V5607-V5603	