

---

## NF-L Mouse mAb, HRP conjugated

Cat. No:	V5508-HRP
Name:	NF-L Mouse mAb, HRP conjugated
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	3A5
Host/Source:	Mouse
Reactivity:	(predicted:&nbsp;Human)
Format:	Liquid
Concentration:	1mg/ml
Buffer:	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Applications:	ELISA
Antigen:	Recombinant Human NF-L Protein