

---

## **phospho-Tau (Thr181) Mouse mAb, Biotin conjugated**

Cat. No:	V5226-Bio
Name:	phospho-Tau (Thr181) Mouse mAb, Biotin conjugated
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	10F4
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	Human
Format:	Liquid
Concentration:	1mg/ml
Buffer:	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Applications:	ELISA
Antigen:	Recombinant Human p-Tau (Thr181) Protein