

---

## Human IgA Mouse mAb

Cat. No:	V20108
Name:	Human IgA Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	9F7
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS(pH7.4)
Applications:	ICA
Antigen:	Purified native protein ( <a href="#">bs-0474P</a> )
Pairing Information:	
V20105-V20108	