
RBP4 Mouse mAb, Biotin conjugated

Cat. No:	V6202-Bio
Name:	RBP4 Mouse mAb, Biotin conjugated
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	7D3
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	(predicted: Human)
Format:	Liquid
Buffer:	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Applications:	ELISA, ICA
Antigen:	Recombinant Human RBP4 Protein