

Recombinant human PDGF-AB protein (Active, CHO)

Catalog Number: bs-48072P

Species: Human

AA Seq: Ser87-Thr211 & Ser82-Thr190

Predicted MW: 14.3 & 12.3

Tags: Tag free

Activity: Yes

Endotoxin: ≤ 10 EU/mg

Purity: $\geq 95\%$ as determined by SDS-PAGE.


Purification: AC

Form: Lyophilized

Storage: Lyophilized from a 0.22 μ m-filtered solution containing PBS, 5% Mannitol and 0.01% Tween 80, pH7.4

Use a manual defrost freezer and avoid repeated freeze-thaw cycles. 36 months at -20°C to -80°C in lyophilized state. 6 months at -20°C to -80°C under sterile conditions after reconstitution. 7-10 days at 2°C to 8°C under sterile conditions after reconstitution.

VALIDATION IMAGES



Measured in a cell proliferation assay using NIH3T3 cells. The ED50 is ≤ 1.0 ng/ml.