
Recombinant human IL36 alpha protein, N-His

Catalog Number: bs-42557P

Concentration: >1mg/ml

Species: Human

AA Seq: 6-158/158

Predicted MW: 19.5 kDa

Tags: N-His

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE


Purification: AC

Form: Liquid

Storage: 20mM Tris-HCl (pH=8.0) with 150mM NaCl

Background: Expressed in immune system and fetal brain, but not in other tissues tested or in multiple hematopoietic cell lines.

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.