
Recombinant Human Tau (acetyl K281), C-His

Catalog Number: bs-47606P

Species: Human

AA Seq: 1-441/441

Predicted MW: 48 kDa

Tags: C-His

Purity: > 85% as determined by SDS-PAGE

Form: Liquid

Storage: 150 mM PBS

Stored at -70°C or -20°C for one year. Avoid repeated freeze/thaw cycles. Stored at 2 to 8 °C for 2 days after thawing.

Background: Tau proteins are important Promotes microtubule assembly and stability, and might be involved in the establishment and maintenance of neuronal polarity. The C-terminus binds axonal microtubules while the N-terminus binds neural plasma membrane components, suggesting that tau functions as a linker protein between both. Axonal polarity is predetermined by tau localization (in the neuronal cell) in the domain of the cell body defined by the centrosome. The short isoforms allow plasticity of the cytoskeleton whereas the longer isoforms may preferentially play a role in its stabilization. Tau proteins subcellular located in the axons of neurons, in the cytosol and in association with plasma membrane components. It expressed in neurons. PNS-tau is expressed in the peripheral nervous system while the others are expressed in the central nervous system.