

**bsm-4485M**

**[ Primary Antibody ]**

## **FEN Mouse mAb**



[www.bioss.com.cn](http://www.bioss.com.cn)

[sales@bioss.com.cn](mailto:sales@bioss.com.cn)

[techsupport@bioss.com.cn](mailto:techsupport@bioss.com.cn)

400-901-9800

### — DATASHEET —

<b>Host:</b> Mouse	<b>Isotype:</b> IgG	<b>Applications:</b> ELISA (1:5000-10000)
<b>Clonality:</b> Monoclonal	<b>CloneNo.:</b> 2F6	
<b>Target:</b> FEN		
<b>Purification:</b> affinity purified by Protein G		