
Human IgE Mouse mAb, FITC conjugated

Cat. No:	V5404-FITC
Name:	Human IgE Mouse mAb, FITC conjugated
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	3D3
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	(predicted: Human)
Format:	Liquid
Concentration:	1mg/ml
Buffer:	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Antigen:	Purified native protein