

bsm-4486M

[Primary Antibody]

COT Mouse mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Mouse

Clonality: Monoclonal

Target: COT

Purification: affinity purified by Protein G

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4).
Shipped at 4°C. Store at -20°C for one year. Avoid repeated
freeze/thaw cycles.

Isotype: IgG

CloneNo.: 7G1

Applications: ELISA (1:5000-10000)

Reactivity: (predicted: COT, Species
independent)