
Human Thyroglobulin Mouse mAb

Cat. No:	V9001
Name:	Human Thyroglobulin Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	1B5
Host/Source:	Mouse
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Applications:	ELISA
Antigen:	Natural human thyroglobulin
Pairing Information:	
V9001-V9002	