

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

Human Thyroglobulin Mouse mAb

Cat. No: V9002

Name: Human Thyroglobulin Mouse mAb

Storage: Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Clone type: Monoclonal

Clone No: 2C1
Host/Source: Mouse
Reactivity: Human

Purity: >90% as determined by SDS-PAGE

Format: Liquid

Buffer: 0.1M PBS(pH=7.4)

Applications: ELISA

Antigen: Natural human throglobulin

Pairing Information:

V9001-V9002