

---

## Aflatoxin B1 Mouse mAb

Cat. No:	V10602
Name:	Aflatoxin B1 Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	4B4
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS (pH7.4)
Applications:	ELISA, ICA
Antigen:	Aflatoxin B1-KLH