
PTH Mouse mAb

Cat. No:	V9101
Name:	PTH Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	4C6
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS(pH7.4)
Applications:	ELISA
Antigen:	Recombinant Human PTH Protein (bs-41496P)
Pairing Information:	
V9101-V9102	