
Human IgE Mouse mAb

Cat. No:	V5408
Name:	Human IgE Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles
Clone type:	Monoclonal
Clone No:	2A8
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Buffer:	PBS (pH7.4)
Applications:	ICA
Antigen:	Purified native protein (bs-2226P)
Pairing Information:	
V5408-V5402	