
Recombinant rabbit CD4 protein, C-His (HEK293)

Catalog Number: bs-43757P

Concentration: >1mg/ml

Species: Rabbit

AA Seq: 26-396/459

Predicted MW: 48 kDa

Tags: C-His

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Storage: Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: This gene encodes a membrane glycoprotein of T lymphocytes that interacts with major histocompatibility complex class II antigens and is also a receptor for the human immunodeficiency virus. This gene is expressed not only in T lymphocytes, but also in B cells, macrophages, and granulocytes. It is also expressed in specific regions of the brain. The protein functions to initiate or augment the early phase of T-cell activation, and may function as an important mediator of indirect neuronal damage in infectious and immune-mediated diseases of the central nervous system. Multiple alternatively spliced transcript variants encoding different isoforms have been identified in this gene. [provided by RefSeq, Aug 2010].

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.