
Recombinant human CD54/ICAM1 protein, C-His (HEK293)

Catalog Number: bs-43764P

Concentration: >0.5mg/ml

Species: Human

AA Seq: 28-480/532

Predicted MW: 51.4 kDa

Tags: C-His

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Liquid

Storage: PBS (pH=7.4)

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: This gene encodes a cell surface glycoprotein which is typically expressed on endothelial cells and cells of the immune system. It binds to integrins of type CD11a / CD18, or CD11b / CD18 and is also exploited by Rhinovirus as a receptor. [provided by RefSeq, Jul 2008].

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.