

---

## Human CEA Mouse mAb, Biotin conjugated

|               |                                                                                 |
|---------------|---------------------------------------------------------------------------------|
| Cat. No:      | V1905-Bio                                                                       |
| Name:         | Human CEA Mouse mAb, Biotin conjugated                                          |
| Storage:      | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| Clone type:   | Monoclonal                                                                      |
| Clone No:     | 7A3                                                                             |
| Host/Source:  | Mouse                                                                           |
| Subtype:      | IgG1                                                                            |
| Reactivity:   | (predicted:&nbsp;Human)                                                         |
| Format:       | Liquid                                                                          |
| Buffer:       | 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.               |
| Applications: | ELISA                                                                           |
| Antigen:      | Recombinant Human CEACAM5 Protein                                               |