
Recombinant human IL12 protein (HEK293)

Catalog Number: bs-43763P

Concentration: >1mg/ml

Species: Human

AA Seq: 23-328/328 & 23-219/219

Predicted MW: 58.9 kDa

Tags: Tag Free

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE


Purification: AC

Form: Liquid

Storage: PBS (pH=7.4)

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.