
Recombinant Chikungunya virus protein (CHIKV) E1 protein, C-His (HEK293)

Catalog Number: bs-43793P

Concentration: >1mg/ml

Species: Chikungunya virus

AA Seq: 810-1224/1248

Predicted MW: 47

Tags: C-His

Endotoxin: Not analyzed

Purity: $\geq 90\%$ as determined by SDS-PAGE.

Purification: AC

Form: Liquid

Storage: PBS (pH=7.4)

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles