

Recombinant rat CD64/FcRI protein, C-His (HEK293)

Catalog Number: bs-43787P

Concentration: >0.5mg/ml

Species: Rat

AA Seq: 16-285/353

Predicted MW: 32.1

Tags: C-His

Endotoxin: Not analyzed

Purity: $\geq 90\%$ as determined by SDS-PAGE.

Purification: AC

Form: Liquid

Storage: PBS (pH=7.4)

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles

Background: This gene encodes a protein that plays an important role in the immune response. This protein is a high-affinity Fc-gamma receptor. The gene is one of three related gene family members located on chromosome 1. [provided by RefSeq, Jul 2008]

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.