

Recombinant human ALPPL2 protein, C-His (HEK293)

Catalog Number: bs-43778P

Concentration: >1mg/ml

Species: Human

AA Seq: 20-503/532

Predicted MW: 54.5

Tags: C-His

Endotoxin: Not analyzed

Purity: $\geq 90\%$ as determined by SDS-PAGE.

Purification: AC

Form: Liquid

Storage: PBS (pH=7.4)

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles

Background: There are at least four distinct but related alkaline phosphatases: intestinal, placental, placental-like, and liver/bone/kidney (tissue non-specific). The intestinal alkaline phosphatase gene encodes a digestive brush-border enzyme. This enzyme is upregulated during small intestinal epithelial cell differentiation. [provided by RefSeq, Jul 2008]

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.