
Recombinant EEEV E2 protein, C-His (HEK293)

Catalog Number: bs-43791P

AA Seq: 323-660/1239

Predicted MW: 39.9

Tags: C-His

Endotoxin: Not analyzed

Purity: $\geq 90\%$ as determined by SDS-PAGE.


Purification: AC

Form: Liquid

Storage: PBS (pH=7.4)

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.