

bs-34296G-BF488

[Secondary Antibodies]

AffiniPure Goat Anti-Mouse IgG, Light Chain Specific, BF488 conjugated

BioSS
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

DATASHEET

Host: Goat

Isotype: IgG

Clonality: Polyclonal

Target: AffiniPure Goat Anti-Mouse IgG, Light Chain Specific

Purification: affinity purified by Protein G

Storage: 10 mM TBS (pH=7.4) with 1% BSA, 0.02% Proclin300 and 50% glycerol.
Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Background: The Fc region is the constant region on an immunoglobulin molecule, the area that is exactly the same on all antibodies. The Fc region is found on the heavy chains and is not involved in binding antigens.

Applications: IF

**Flow-Cyt
ICC/IF**

Excitation Spectrum: 488nm

Emission spectrum: 519nm

Reactivity: Mouse