

bs-34296G-BF680

[Secondary Antibodies]

AffiniPure Goat Anti-Mouse IgG, Light Chain Specific, BF680 conjugated

BioSS
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Goat

Isotype: IgG

Clonality: Polyclonal

Target: AffiniPure Goat Anti-Mouse IgG, Light Chain Specific

Purification: affinity purified by Protein G

Storage: 10 mM TBS (pH=7.4) with 1% BSA, 0.02% Proclin300 and 50% glycerol.
hipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Background: The Fc region is the constant region on an immunoglobulin molecule, the area that is exactly the same on all antibodies. The Fc region is found on the heavy chains and is not involved in binding antigens.

Applications: IF

**Flow-Cyt
ICC/IF**

Excitation Spectrum: 680nm

Emission spectrum: 702nm

Reactivity: Mouse