
Human HAVCR1/KIM-1 Mouse mAb

Cat. No:	V8002
Name:	Human HAVCR1/KIM-1 Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	12A11
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Concentration:	1mg/ml
Buffer:	0.01M PBS(pH7.4)
Applications:	ELISA
Antigen:	Recombinant Human HAVCR1 Protein (bs-43178P)
Pairing Information:	
V8001-V8002	