
Human Vasopressin-neurophysin II-copeptin Mouse mAb

Cat. No:	V9201
Name:	Human Vasopressin-neurophysin II-copeptin Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	8G11
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Concentration:	1mg/ml
Buffer:	0.01M PBS(pH7.4)
Applications:	ELISA
Antigen:	Recombinant Human Vasopressin-neurophysin 2-copeptin Protein (bs-43179P)
Pairing Information:	V9201-V9202&V9203