
Human Vitamin D Binding Protein Mouse mAb

Cat. No:	V6902
Name:	Human Vitamin D Binding Protein Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	12C9
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Concentration:	1mg/ml
Buffer:	0.01M PBS(pH7.4)
Applications:	ELISA, ICA
Antigen:	Recombinant Human Vitamin D Binding Protein (bs-43119P)
Pairing Information:	V6901&V6902-V6903