

Recombinant mouse CD28 protein, C-His (HEK293)

Catalog Number: bs-43795P

Concentration: >1mg/ml

Species: Mouse

AA Seq: 21-149/218

Predicted MW: 16.8

Tags: C-His

Endotoxin: Not analyzed

Purity: $\geq 90\%$ as determined by SDS-PAGE.

Purification: AC


Form: Liquid

Storage: PBS (pH=7.4)

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles

Background: The protein encoded by this gene is essential for T-cell proliferation and survival, cytokine production, and T-helper type-2 development. Several alternatively spliced transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Jul 2011]

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.