
Recombinant human FTMT protein, N-His

Catalog Number: bs-42601P

Concentration: >1mg/ml

Species: Human

AA Seq: 50-242/242

Predicted MW: 39

Tags: N-His

Endotoxin: Not analyzed

Purity: $\geq 90\%$ as determined by SDS-PAGE.

Purification: AC

Form: Liquid

Storage: 20MMTris-HCl (pH=8.0) with 4M Urea, 150mM NaCl

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.