
Recombinant human ANX8 protein, N-His

Catalog Number: bs-42588P

Species: Human

Tags: N-His

Purity: $\geq 90\%$ as determined by SDS-PAGE.


Purification: AC

Form: Liquid

Storage: 20MMTris-HCl (pH=8.0) with 150mM NaCl

Background: Annexin A8 like 2 is a member of the annexin family of evolutionarily conserved Ca^{2+} and phospholipid binding proteins. The protein may function as an anticoagulant that indirectly inhibits the thromboplastin-specific complex, and overexpression of this gene has been associated with acute myelocytic leukemia. A highly similar duplicated copy of this gene is found in close proximity on the long arm of chromosome 10.

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.