

---

## ARFIP1 Rabbit pAb

Catalog Number: bs-1477R

Target Protein: ARFIP1

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Pig, Cow, Chicken, Dog, Horse)

Predicted MW: 41 kDa

Entrez Gene: 27236

Swiss Prot: P53367

Source: KLH conjugated synthetic peptide derived from human Arfaptin 1: 255-373/373.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: human homolog binds ADP-ribosylation factor (ARF) and inhibits ARF dependent phorbol ester induced secretion of matrix metalloproteinase 9 [RGD]