

bs-10260R**[Primary Antibody]****BCAS4 Rabbit pAb**

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human) Predicted MW.: 23 kDa Subcellular Location: Cytoplasm
Clonality: Polyclonal		
GeneID: 55653	SWISS: Q8TDM0	
Target: BCAS4		
Immunogen: KLH conjugated synthetic peptide derived from human BCAS4: 61-160/211.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: BCAS4 (Breast Carcinoma Amplified Sequence 4) is a Protein Coding gene. Diseases associated with BCAS4 include choanal atresia. An important paralog of this gene is BLOC1S4. BCAS4 A chromosomal aberration involving BCAS4 has been found in some breast carcinoma cell lines. Translocation t(17;20)(q23;q13) with BCAS3. Belongs to the cappuccino family. 3 isoforms of the human protein are produced by alternative splicing.		