

bs-18080R

[Primary Antibody]

HSPA2 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Pig, Sheep, Cow, Dog, Horse, Monkey) Predicted MW.: 70 kDa Subcellular Location: Cytoplasm
Clonality: Polyclonal		
GeneID: 3306	SWISS: P54652	
Target: HSPA2		
Immunogen: KLH conjugated synthetic peptide derived from human HSPA2: 451-550/639.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		

— SELECTED CITATIONS —

- **[IF=13.751]** Jianjun Hu. et al. Ultrasound combined with nanobubbles promotes systemic anticancer immunity and augments anti-PD1 efficacy. J Immunother Cancer. 2022 Mar;10(3):e003408 WB ;Mouse. 35236741
- **[IF=2.942]** Huang Y et al. Ubiquitome analysis reveals the involvement of lysine ubiquitination in the spermatogenesis process of adult buffalo (Bubalus bubalis) testis. Biosci Rep. 2020 Jun 26;40(6):BSR20193537. WB ;swamp buffalo (Bubalus bubalis). 32469046
- **[IF=2.823]** Peng - fei Zhang. et al. Integrated analysis of phosphoproteome and ubiquitylome in epididymal sperm of buffalo (Bubalus bubalis). Mol Reprod Dev. 2021 Jan;88(1):15-33 WB,IP ;Bovine. 33140506